

QX WORLD®

QUEX PET PADS USER GUIDE:

ADVANCED COMFORT AND WELLBEING FOR YOUR PETS

pet:pad

Introduction

Transform your pet's wellness routine with the new QUEX Pet Pads. Designed with cutting-edge QUEX World technology, these pads offer superior conductive performance to enhance your pet's wellbeing and comfort. Suitable for both small and large animals, QUEX Pet Pads provide enhanced energetic coverage through four large, highly conductive carbon-impregnated rubber electrodes.

Key Features:

- Enhanced Effectiveness: Four large electrodes provide extensive energetic coverage for optimal wellbeing.
- High Conductivity & Durability: Made with carbon-impregnated rubber, the pads are built to last, delivering consistent performance with every use.
- Versatile Connectivity: Seamlessly connects to the limb port at the front of the Quex S or back of the Quex ED with a 7.2-meters cable (1.2 meters permanently attached + 6 meters extension).
- Portable Convenience: The pads easily fit into a carry case for on-the-go sessions or can be placed on your pet's favorite resting spot.

Sizes:

Available in two sizes to accommodate both small and large animals, ensuring that every pet can benefit from advanced biofeedback support.

Small Pet Pad:

Pad Size: 65 x 45 cm
Individual Electrode Size: 30 x 20 cm
Ideal For: Cats, small dogs (up to 20 kg), hamsters, guinea pigs

Large Pet Pad:

Pad Size: 84 x 84 cm
Individual Electrode Size: 40 x 40 cm
Ideal For: Large dogs (above 20 kg), horses, cattle (can be placed on the back)

Connection Details:

Cable Length: 1.2 meters with an additional 6-meters extension for maximum flexibility.
Compatible devices: QUEX ED, QUEX S, SCIO, EDUCTOR (additional connector is required for SCIO).

Preparation for Use

The rubber material is protected by a thin layer of cling, which should be carefully removed before use. The film should peel off easily. The wiring, which is the sensitive part of the equipment, is housed underneath the Velcro. To prevent damage, avoid handling the Velcro strips on the sides.

Safety Precautions

Ensure that the plug connections are kept secure and dry, especially when the device is used around small children and animals.

Care Instructions

To clean the product, use a soft, damp cloth with a small amount of gentle soap solution if needed. For soiled areas, spray 100% AM Alcohol spray on a cloth and use it to wipe the pad down.

Avoid using harsh chemicals, and not submerge the accessory in water or any liquid. If accidental submersion occurs, allow the pad to dry out completely and test it before further use.

Give your pets the gift of advanced care and comfort with QUEX Pet Pads. Enhance their recovery, performance, and overall health!



QX WORLD®

QX WORLD LTD.

Hungary
1095 Budapest, Tinodi Street 1-3, building A, floor IV, door 93
E-mail: info@qxworld.eu

www.qxworld.eu

QX World Ltd. products are not medical devices, these are not designed or intended for use in the diagnosis, treatment, mitigation, or prevention of any disease, medical condition, or disorder. The Product should not be used as a substitute for professional medical advice, diagnosis, or treatment. Users should always consult their healthcare provider before making lifestyle changes or in cases of health concerns. The device was tested according to the stricter medical device requirements. The device is under continuous development with the aim of becoming a certified diagnostic medical device.

This document is protected by copyright. You must not copy, reproduce, distribute, perform, broadcast or otherwise make available for others, exhibit, adapt or otherwise change this document unless you receive a prior written license from QX World Kft. Any unauthorized use is strictly prohibited and QX World Kft. reserves all rights to enforce its copyright before the competent courts or authorities.